

Product Information

FreezeMe One, Cryopreservation Medium with Fetal Bovine Serum
 Cat. No. FM1-F (50 ml)

General Information

FreezeMe One is a cryoprotective medium that contains 10 % DMSO and 20 % FBS. The serum used is extensively tested to protect cells during cell preservation.

Features:

- Cryopreservation of a wide range of cell types with high viability
- Ready-to-use solution

Product Specifications

Appearance	Clear red liquid
Storage and shelf life	Store at $\leq -15^{\circ}\text{C}$. FreezeMe One is a light sensitive solution. It should be protected from light during storage.
Shipping conditions	Frozen (Dry Ice)
Thawing	$+37^{\circ}\text{C}$ water bath or overnight at $+2^{\circ}\text{C}$ to $+8^{\circ}\text{C}$. Swirl gently to homogenize.

Quality Control

Freezing Protocols:

Before cryopreservation cells should be checked for contamination. FreezeMe One can be used with any standard freezing protocol.

Cryopreservation of Suspension Cultures:

1. Count the number of viable cells to be cryopreserved. Cells should be in mid-log phase of growth. Centrifuge the cells for 5 min to pellet cells (200 to 400 g). Remove the supernatant down to the smallest volume without disturbing the cells.
2. Resuspend cells in pre-cooled ($+4^{\circ}\text{C}$ to $+8^{\circ}\text{C}$) FreezeMe One to a concentration of 5×10^6 to 1×10^7 cells/ml.
3. Aliquot into cryogenic storage vials. Place vials at $+4^{\circ}\text{C}$ and start the freezing procedure within 5 min. Cells are frozen slowly at $+1^{\circ}\text{C}/\text{min}$ (by programmable coolers or by placing vials in an insulated box in a -70°C to -90°C freezer).
4. Then transfer storage vials to liquid nitrogen storage.

Cryopreservation of Adherent Cultures:

1. Detach cells from the substrate with a gentle dissociating agent. Especially with sensitive cells use Accutase* (Cat. No. ACC-1B) to avoid cell damage. Inactivate dissociating agent if necessary.
2. Resuspend the detached cells in complete growth medium and establish the viable cell count.
3. Centrifuge for 5 min to pellet cells (200 to 400 g). Remove the supernatant down to the smallest volume without disturbing the cells.
4. Resuspend cells in pre-cooled ($+4^{\circ}\text{C}$ to $+8^{\circ}\text{C}$) FreezeMe One to a concentration of 5×10^6 to 10^7 cells/ml.

Product Information

FreezeMe One, Cryopreservation Medium with Fetal Bovine Serum

Cat. No. FM1-F (50 ml)

5. Aliquot into cryogenic storage vials. Place vials at +4°C and start the freezing procedure within 5 min. Cells are frozen slowly at +1°C/min (by programmable coolers or by placing vials in an insulated box in a -70°C to -90°C freezer).
6. Then transfer storage vials to liquid nitrogen storage.

Thawing of Cryopreserved Cells:

Cryopreserved cells can be thawed by the following procedures:

Centrifugation:

1. Remove cells from storage and thaw quickly in a +37°C water bath. Capricorn Scientific recommends eye protection by using approved safety goggles. We also suggest the use of safety gloves to protect uncovered skin.
2. Place 1 to 2 ml of thawed cells in ~25 ml of complete growth medium. Mix cell suspension gently.
3. Centrifuge the cells at ~80 g for 2 to 3 min.
4. Check clarity of the supernatant and visibility of a consolidated cell pellet. Discard supernatant without disturbing the cells.
5. Gently resuspend the cells in complete growth medium and perform a viable cell count.
6. Plate the cells. Cell inoculum should be at least 3×10^5 viable cells/ml.

Direct Plating:

1. Remove cells from storage and thaw quickly in a +37°C water bath. Capricorn Scientific recommends eye protection by using approved safety goggles. We also suggest the use of safety gloves to protect uncovered skin.
2. Plate cells directly with complete growth medium. Use 10 to 20 ml of complete medium per 1 ml of frozen cells. Cell inoculum should be at least 3×10^5 cells/ml.
3. Culture cells for 12 to 24 h. Replace medium with fresh complete growth medium to remove cryopreservative residues.

We recommend thawing procedure 1, especially when handling sensitive cells.

Precautions and Disclaimer

This product is for research use and further manufacturing only.

Help Needed?

If you have any further questions regarding this product, please do not hesitate to contact our cell culture experts by email (techservice@capricorn-scientific.com) or phone (+49 6424 944640).

*Accutase is a registered trademark of Innovative Cell Technologies, Inc.