



# qPCR TROUBLESHOOTING GUIDE



appleton<sup>®</sup>



| Problem  | Possible Cause   | Recommendation   |
|--|--|--|
| <b>No fluorescent signal trace and no product on agarose gel</b> | Activation time too short  | Ensure qPCR mix is activated for a minimum of 2min at 95°C before cycling.   |
|  | Error in protocol setup  | Verify that correct reagent concentrations, volumes, dilutions and storage conditions have been used.  |
|  | Suboptimal primer/probe design   | Use primer design software or validated primers/probes. Test assay on a control template.  |
|  | Incorrect concentration of primers/probe   | Use primer concentrations between 200nM and 1µM. To optimize probe concentration, test the probe at several levels from 50 to 250 nM final concentrations in PCR with optimized levels of primers; generally, the probe should be at least 2-fold lower than the primer concentration. |
|  | Template degraded  | Re-isolate your template from the sample material or use freshly prepared template dilution. Verify the integrity of RNA using agarose gel electrophoresis.  |
|  | Primers/probe degraded   | Use newly synthesized primers and/or probe   |
|  | Template contaminated with qPCR inhibitors   | Perform a dilution series of the PCR template to determine whether the effect of the inhibitory agent can be reduced.  |
|  |  | Re-purify template and re-suspend it in DEPC-treated water.  |
|  | Template concentration too low   | Increase concentration used - see protocol for ideal concentrations depending on type of template.   |
| Cycling conditions not optimal                                   | Increase extension/annealing time, increase cycle number, reduce annealing temperature |  |





| Problem  | Possible Cause                                 | Recommendation   |
|--|--|--|
| <p><b>PCR product present on agarose gel but no fluorescent signal trace</b></p> | <p>Error in instrument setup</p>               | <p>Check that the acquisition settings are correct during cycling</p>  |
| <p><b>Non-specific amplification and/or primer-dimers</b></p>                    | <p>Inefficient reverse transcription</p>       | <p>Extend reverse transcription time up to 30min and/or increase the temperature up to 55°C</p>  |
|  | <p>Suboptimal primer/probe design</p>          | <p>Redesign primers and/or probe using appropriate software, or use validated primers/probes. Perform melt-curve analysis to check if primer-dimers are present / absent.</p>  |
|  | <p>Suboptimal primer/probe concentration</p>   | <p>Use primer concentration between 100nM and 1µM. To optimize probe concentration, test the probe at several levels from 50 to 250 nM final concentrations in PCR with optimized levels of primers; generally, the probe should be at least 2-fold lower than the primer concentration.</p> |
|  | <p>Annealing/extension temperature too low</p> | <p>Increase annealing/extension temperature up to 65°C or until primer-dimer/non-specific amplification products disappear.</p>  |
|  | <p>Suboptimal template concentration</p>       | <p>Vary template concentration until non-specific products disappear</p>   |
|  | <p>Extension time too long</p>                 | <p>Reduce extension time to determine whether non-specific products are reduced.</p>   |



| Problem   | Possible Cause   | Recommendation  |
|---|--|---|
| Late amplification trace  | Inefficient reverse transcription  | Extend reverse transcription time up to 30min and/or increase the temperature up to 55°C.   |
|   | Activation time too short  | Ensure qPCR mix is activated for 2-3mins at 95°C before cycling, depending on type of template used - see protocol.   |
|   | Annealing temperature too high   | Optimise annealing temperature in steps of 2°C using a thermal gradient.  |
|   | Extension time too short   | Double extension time to determine whether the cycle threshold (CT) is affected.  |
|   | Template concentration too low   | Increase concentration depending on template type - see protocol.   |
|   | Degraded or sheared nucleic acid template  | Re-isolate your template from the sample material or use freshly prepared template dilution. Confirm the integrity of RNA using agarose gel electrophoresis.  |
|   |  | Do not store diluted template in water or at low concentrations.  |
|   | Suboptimal primer/probe design   | Evaluate primer sequences for complementarity and secondary structure. Redesign primers/probe using appropriate software, or use validated primers.   |
|   | Suboptimal primer/probe concentration  | Use primer concentration between 100nM and 1µM. To optimize probe concentration, test the probe at several levels from 50-250nM final concentrations in PCR with optimized levels of primers; generally, the probe should be at least 2-fold lower than the primer concentration. |
|   | Assay Design - PCR target may be too long or may contain too much secondary structure leading to inefficient PCR | qPCR is most efficient when the PCR product is smaller in length. Select a target sequence with between 80 - 200bp for optimal results.   |
| Follow these guidelines to avoid secondary structure: Use an annealing temperature above the melting temperature (T <sub>m</sub> ) of any template secondary structure Avoid templates with long (>4) repeats of single bases Maintain a GC content of 50–60% Analyse secondary structure using the DNA mfold server created by Dr Michael Zuker ( <a href="http://www.bioinfo.rpi.edu/applications/mfold">http://www.bioinfo.rpi.edu/applications/mfold</a> ) or an equivalent primer-design program |  |   |



| Problem   | Possible Cause                      | Recommendation   |
|---|-------------------------------------|--|
| <b>High background signal or gradual increase in background fluorescence due to inefficient fluorophore quenching</b> | Partially degraded / unstable probe | Do not store diluted probes in water or at low concentrations.   |
|   |                                     | Do not store hybridisation probes at temperatures above $-20^{\circ}\text{C}$ (probes are prone to degradation at $> -20^{\circ}\text{C}$ , as bonds between the oligonucleotide probe and the conjugated fluorophore become labile).  |
|   |                                     | Check with the oligonucleotide manufacturer for the recommended dilution and storage conditions for the hybridisation probe.   |
|   |                                     | Aliquot probe stock into small volumes and thaw each aliquot only once (minimises multiple freeze-thaw cycles which degrade probes)  |
| <b>Signal in negative control</b>   | Poor laboratory technique           | To minimise the possibility of contamination of PCR components by PCR product or other template, designate a work area exclusively for PCR assay setup.  |
|   |                                     | Use a solution of 10% bleach instead of ethanol to prepare the workstation area for PCR assay setup. Ethanol will only induce precipitation of DNA in your work area, while the 10% bleach solution will hydrolyse, as well as dissolve, any residual DNA.   |
|   |                                     | In general, follow these practices to minimize the risk of sample contamination: <ul style="list-style-type: none"><li>- Wear gloves</li><li>- Use screwcap tubes</li><li>- Use aerosol-resistant filter tips</li><li>- Use calibrated pipets dedicated to PCR</li><li>- Use PCR-grade water and use it only for PCR</li><li>- Use a no-template control to verify absence of contamination</li><li>- Prepare a master mix with sufficient volume to prepare all replicate samples</li></ul> |





| <b>Problem</b>  | <b>Possible Cause</b>  | <b>Recommendation</b>  |
|---|--|--|
| <b>Amplification of genomic DNA in no RT-control</b>          | Genomic DNA contamination of RNA used in RT-qPCR   | Treat samples with purified RNase-free DNase before reverse transcription.   |
|   |  | It may also be helpful to design primers at splice junctions to avoid genomic DNA amplification.   |
| <b>Poor CT value reproducibility across replicate samples</b> | Poor laboratory technique or imprecise pipetting   | In general, follow these practices to improve replicate reproducibility: <ul style="list-style-type: none"><li>- Use aerosol-resistant filter tips</li><li>- Use calibrated pipets dedicated to PCR</li><li>- Use PCR-grade water and use it only for PCR</li><li>- Prepare a master mix with sufficient volume to prepare all replicate samples</li><li>- Add template to master mix before pipetting into reaction vessels</li><li>- Avoid pipetting less than 5µl</li></ul> |
|   | Inhibitors of PCR are carried over from sample preparation   | Perform a dilution series of the PCR template to determine whether the effect of the inhibitory agent can be reduced.<br>Take extra care with the nucleic acid extraction steps to minimize carryover of PCR inhibitors like phenol, detergents and proteases.   |
|   | Some primers are particularly sensitive to thermal cycling conditions, leading to poor reproducibility in amplification reactions. | Determine how your primers behave at different annealing temperatures by doing a gradient PCR.   |
|   | Air bubbles in reaction mix  | Centrifuge reaction samples/plate prior to running on qPCR thermal cycler.   |

