

### TIGER STRIPE

18 01 04,  
20 01 99  
& 18 02 03



Available in  
1L, 2L,  
5L, 6L  
& 30L

Offensive  
Hygiene Waste

Deep Landfill/  
Energy From  
Waste

### ORANGE

18 02 02\*  
& 18 01 03\*



Available in  
1L, 2L,  
5L, 6L  
Loop & 30L

Infectious  
Waste

For Alternative  
Treatment

### YELLOW

18 01 03\*  
& 18 01 09



Available in  
1L, 2L,  
5L, 6L  
Loop & 30L

Infectious  
Waste

For Incineration  
Only

### RED

18 01 02,  
18 01 03\*  
& 18 02 02\*



Available in  
1L, 2L,  
5L & 30L

Anatomical  
Waste

For Incineration  
Only

### PURPLE

18 01 03\*  
& 18 02 08



Available in  
1L, 5L, 6L  
& 30L

Cytotoxic Or  
Cytostatic  
Contaminated  
Waste

For Incineration  
Only

### BLUE

18 01 09  
& 18 02 08



Available in  
1L, 5L  
& 30L

Medicinal  
Waste

For Incineration  
Only

# PRODUCTS & ACCESSORIES

Reduce your carbon footprint by up to **90%** compared to plastic waste bins



Helping you save money and the environment

## 96% PAPER BASED

*non-sharp clinical waste containers -  
The alternative to plastic for  
all your non-sharp waste!*



**Bio-bin®** Loop, 6 Litre and 30 Litre **Bio-bin®**

also available on special order as a green and white bin without any waste labelling

[www.bio-bin.co.uk](http://www.bio-bin.co.uk)

TEL: +44 (0) 114 327 7573

EMAIL: [info@bio-bin.co.uk](mailto:info@bio-bin.co.uk)



Bio-bin® is a UN 3291 approved waste container for non-sharp medical and clinical waste

## SIZE RANGE



**Bio-bin®**  
1 Ltr  
1 Litre  
Max fill weight 0.5kg  
Case Size: 100 bins



**Bio-bin®**  
2 Ltr  
2 Litre  
Max fill weight 1kg  
Case Size: 100 bins

### Carbon Footprint



**Bio-bin®**  
Loop  
Loop  
Max fill weight 2kg  
Case Size: 100 bins



**Bio-bin®**  
5 Ltr  
5 Litre  
Max fill weight 2kg  
Case Size: 100 bins  
0.18KG CO2



**Bio-bin®**  
6 Ltr  
6 Litre  
Max fill weight 1kg  
Case Size: 40 bins  
0.19KG CO2



**Bio-bin®**  
30 Ltr  
30 Litre  
Max fill weight 8kg  
Case Size: 10 bins  
1.22KG CO2

ALL OUTER CASES MADE FROM 100% RECYCLED MATERIAL

## Bio-bin® PLACENTA-PACK

Each placenta pack contains 100 plastic zip-loc bags and 50 5-Litre red Bio-bins to collect your placenta waste.



### SIZE AVAILABILITY FULL CASES

	1LTR	2LTR	LOOP	5LTR	6LTR	30LTR
Qty per case	100	100	100	100	40	10
YELLOW	YES	YES	YES	YES	YES	YES
RED	YES	YES		YES		YES
ORANGE	YES	YES	YES	YES	YES	YES
BLUE	YES			YES		YES
TIGER STRIPE	YES	YES		YES	YES	YES
PURPLE	YES			YES	YES	YES
PLAIN			YES		YES	YES

### SIZE AVAILABILITY MINIPACKS

	1LTR	2LTR	5LTR
Qty per case	20	20	20
YELLOW	YES	YES	YES
RED	YES	YES	YES
ORANGE	YES	YES	YES
BLUE	YES		YES
TIGER STRIPE	YES	YES	YES
PURPLE	YES		YES

## 1L, 2L & 5L Bio-Bin Mini-Packs 20 Bio-Bins per case


Perfect for smaller units, community nursing & labs.



## ACCESSORIES



M Stand for 2 & 5 Litre Bio-bin®



Loop M Stand Bio-bin®



Z Stand for 6 Litre Bio-bin®



U Stand for 6 Litre Bio-bin®




Trolley for 30 Litre Bio-bin®



Wheeled Floor Stand for 30 Litre Bio-bin®



M Stand for 30L Bio-bin®



Floor Stand for 30 Litre Bio-bin®