



**PRODUCT INFORMATION**

<b>Item Number</b>	<b>GL890</b>
<b>Product Colour</b>	<b>Blue</b>
<b>Size Range</b>	<b>XS-XXL</b>
<b>Fingertip Thickness</b>	<b>0.11 mm</b>
<b>Palm Thickness</b>	<b>0.08 mm</b>
<b>Typical Weight</b>	<b>5.00 g</b>
<b>Shelf Life</b>	<b>5 Years</b>

**FEATURES**

**Chemical Resistant:** Achieves type C chemical resistance according to EN ISO 374-1, for chemical protection against at least 1 chemical

**EN ISO 374-5 Protection:** Conforms with the requirements of EN ISO 374-5, for protection against bacteria, fungi and virus

**Sensitivity:** Latex free, reducing the risk of hypersensitivity and allergic reactions to natural rubber proteins

**Grip:** Textured finish to the fingers to enhance grip and handling

**Food Safe:** Suitable for food contact; all food types. Please refer to our online Food Compatibility Chart for full details

**Comfort:** Nitrile technology provides an exceptionally soft and pliable material

**Material Composition:** Manufactured from nitrile butadiene polymer

**Powder Content:** Powder Free, reducing the risk of dust contamination

**Silicone Free:** Silicone residues from gloves can interfere with painting, coating, or adhesive bonding processes, potentially causing defects. These silicone-free gloves help minimise that risk, ensuring clean application surfaces.

**SECTORS**

Automotive • Engineering and Manufacturing • Food Manufacturing and Preparation • Janitorial • Medical, Vet, and Dental



**STANDARDS**

Tested to European standards EN ISO 21420, EN ISO 374-1, EN ISO 374-4 & EN ISO 374-5

AQL 1.5

EU Certified

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425

Compliant with Medical Devices Regulation (EU) 2017/745 Class 1

**Compliant with European Standards EN455:** Medical Gloves for Single Use parts 1,2,3&4

Tested in compliance with Commission Regulation (EU) No 10/2011 for Food contact



MDR - Class I  
PPE CAT III



2777

VIRUS

AQL



Disposable Gloves

Nitrile

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555



**CAUTION**

Contains accelerator chemicals which may cause type IV skin irritation and contact dermatitis in sensitised individuals.

**PACKING**

100 gloves per dispenser  
10 dispensers per case  
50 cases per pallet

**INNER PACKAGING**



**STORAGE INSTRUCTIONS**

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

**DISPOSAL INSTRUCTIONS**

Dispose of the product responsibly in accordance with local regulations.

**OUTER PACKAGING**



Size	Inner Packaging Barcode
XS	5010699528322
Small	5010699528339
Medium	5010699528346
Large	5010699528353
XL	5010699528360
XXL	5010699536570

Size	Outer Packaging Barcode
XS	5010699839336
Small	5010699839343
Medium	5010699839350
Large	5010699839367
XL	5010699839374
XXL	5010699849809

Inner Packaging	
Inner dimensions (LxWxH)	235 x 117 x 78 mm

Outer Packaging	
Outer Dimensions (LxWxH)	418 x 258 x 251mm

**PLEASE NOTE:**

Weight does not form part of the specification and is subject to change. Palm thickness is approximate and based on an agreed manufacturing tolerance. All thicknesses and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555