



- Extensive Range
- 90+ Years of History
- Extended Aftercare Support
- UKAS Laboratory

Sparkfree Freezers

Sparkfree Freezer RLVF1517

The Labcold Sparkfree range has a specially manufactured sparkfree interior with all sources of potential ignition removed from the chamber. Anything that could cause a spark is located outside of the chamber making it an important part of laboratory safety. For extra peace of mind this freezer is backed by a 2 year warranty.

For extra security a built in lock is fitted. The temperature display features high and low audio and visual temperature alarms along with a door open alarm.

Specifications

Product Code	RLVF1517
Style	Freestanding
Gross Capacity	406 litres
UK Warranty	2 years
External Dimensions (mm) HxWxD	1875mm x 600mm x 700mm
Internal Dimensions (mm) HxWxD	1690mm x 440mm x 500mm
Net Weight	97kg
Shelves	6


Specifications (cont.)

Temperature Range	-18°C to -25°C
Factory Set Point	-20°C
Door Style	Solid
Lock Type	Key lock
Temperature Display	Yes
Self Closing Door	Yes
High/Low Temperature Alarm	Yes
Door Open Alarm	Yes
Mains Failure Alarm	No
Access Porthole	Yes
Access Porthole Size	12mm
Manual Defrost	Yes
Power Consumption Kwh/24hrs	3.16
Eco Friendly Refrigerant	Yes
Ambient Operating Temperature	16°C to 32°C

Features

- Sparkfree safety interior
- Door Lock
- Self closing door
- High/Low temperature alarm
- 6 shelves
- Manual defrost
- Access Port
- Temperature display
- Door Open alarm
- 2 Year parts and labour warranty (UK only)

Models

Model	Gross Capacity	Door Style	Lock Type
 Sparkfree Freezer RLVF1517	406 litres	Solid	Key lock