



APPLETON WOODS LTD
EQUIPMENT & CONSUMABLES
FOR LIFE SCIENCE LABS
www.appletonwoods.co.uk

appleton®

For the Life Scientist


NEW Pipette Tips



High quality, eco-friendly, filter and non-filter pipette tips available to pre-order now



Eco-friendly refills
Low retention
Universal fitting
Filter and non-filter
Sterile
Cardboard rack option



Key Features

Eco friendly

- Refill racks for both filter and non-filter tips available
- Recyclable cardboard packaging used for all refills - reduces plastic waste by 70%
- Tip box recycling available


Low Retention

Less fluid retention for maximum accuracy

- No silicon additives so no leaching
- Specially designed resins create hydrophobic internal surface

Filter Tips

- Racked and refills available
- PE filters free from chemical additives
- Aerosol-tight filters



Filter Tips

10uL	XL, sterile, racked	960 tips	ACA619
	XL, sterile, reload decks	960 tips	ACA671
20uL	Low retention, sterile, racked	960 tips	ACA626
	Low retention, sterile, reload decs	960 tips	ACA672
100uL	Low retention, sterile, racked	960 tips	ACA628
	Low retention, sterile, reload decs	960 tips	ACA673
200uL	Low retention, sterile, racked	960 tips	ACA659
	Low retention, sterile, reload decs	960 tips	ACA674
1250uL	Low retention, sterile, racked	576 tips	ACA668
	Low retention, sterile, reload decs	576 tips	ACA675

Standard Tips

10uL	XL, low retention, sterile, racked	960 tips	ACA613
	XL, sterile, reload decks	960 tips	ACA611
200uL	Low retention, sterile, racked	960 tips	ACA623
	Low retention, sterile, reload decks	960 tips	ACA624
1250uL	Low retention, sterile, racked	576 tips	ACA663
	Low retention, sterile, reload decks	576 tips	ACA664

