

## Product Information

Minimal Essential Medium (MEM) with Earle's Salts, with Stable Glutamine  
 Cat. No. MEM-STA (500 ml)

### Product Specifications

Appearance	Clear red orange solution
CO <sub>2</sub> concentration, optimum	5.0 %
Storage and shelf life	Store at +2°C to +8°C protected from light. Once opened, store at 4° C and use within 6-8 weeks.
Shipping conditions	Ambient

### Formulation

Components	Concentration mg/L
<b>Amino Acids:</b>	
L-Alanyl-L-Glutamine	434.00
L-Arginine HCl	126.00
L-Cystine 2 HCl	31.30
L-Histidine HCl H <sub>2</sub> O	42.00
L-Isoleucine	52.00
L-Leucine	52.00
L-Lysine HCl	72.50
L-Methionine	15.00
L-Phenylalanine	32.00
L-Threonine	48.00
L-Tryptophan	10.00
L-Tyrosine 2 Na 2 H <sub>2</sub> O	51.90
L-Valine	46.00
<b>Vitamins:</b>	
Choline chloride	1.00
D-Calcium Pantothenate	1.00
Folic Acid	1.00

Components	Concentration mg/L
myo-Inositol	2.00
Nicotinamide	1.00
Pyridoxal HCl	1.00
Riboflavin	0.10
Thiamine HCl	1.00
<b>Inorganic Salts:</b>	
CaCl <sub>2</sub> 2 H <sub>2</sub> O	265.00
MgSO <sub>4</sub>	97.67
KCl	400.00
NaHCO <sub>3</sub>	2200.00
NaCl	6800.00
NaH <sub>2</sub> PO <sub>4</sub> (anhydrous)	122.00
<b>Other Components:</b>	
D-Glucose	1000.00
Phenol Red Sodium Salt	11.0

### Precautions and Disclaimer

This product is for research and further manufacturing use only. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

### Help Needed?

If you have any further questions regarding this product, please do not hesitate to contact our cell culture experts by email (techservice@capricorn-scientific.com) or phone (+49 6424 944640).