



## UNDER-COUNTER LABORATORY FRIDGE-FREEZER COMBI SOLID DOOR 119L WHITE



- ✓ **Operating Temperature Range of +2° C to +10° C for Fridge and -18° C to -25° C for Freezer**  
Conforms to MHRA guidelines
- ✓ **NEW Advanced Controller**  
Records temperatures accurate to 0.1° C
- ✓ **NEW Unique Anti Microbial Handle**  
Helps prevent the growth of bacteria and minimise the potential for cross contamination
- ✓ **Fitted Lock as Standard**  
Maintains security and ensures the unit can only be accessed by approved personnel
- ✓ **Fitted Rollers**  
For ease of movement
- ✓ **Wire Shelves**  
Enables optimum air flow inside unit



[www.facebook.com/Lec-Medical](https://www.facebook.com/Lec-Medical)  
[www.twitter.com/Lec\\_Medical](https://www.twitter.com/Lec_Medical)

SKU CODE:	PRODUCT CODE:	DESCRIPTION:
444440972	LSC119UK	Under-counter Laboratory Fridge-freezer Combi Solid Door 119L White



## FULL SPECIFICATION – LSC119UK UNDER-COUNTER LABORATORY FRIDGE-FREEZER COMBI SOLID DOOR 119L WHITE



<b>Description</b>	Under-counter
<b>SKU Code</b>	444440972
<b>EAN Code</b>	5052263009721
<b>GDPA Code</b>	LSC119UK
<b>Gross Capacity (L)</b>	119 (Fridge: 106, Freezer: 13)
<b>Net Capacity (L)</b>	117
<b>Unit External Width (mm)</b>	540
<b>Unit External Depth (mm)</b>	600
<b>Unit External Height (mm)</b>	860
<b>Product Weight (kg)</b>	34
<b>Packed Width (mm)</b>	602
<b>Packed Depth (mm)</b>	639
<b>Packed Height (mm)</b>	877
<b>Packed Weight (kg)</b>	36
<b>Country of Origin</b>	UK
<b>Internal Dimensions (mm) (WxDxH)</b>	Fridge: 450 x 450 x 525, Freezer: 400 x 295 x 120 (step: 155d x 225h)
<b>External Finish</b>	White
<b>Internal Finish</b>	White
<b>Type of Plastic</b>	SRL800 (HIPS/ABS Mix)
<b>Temperature Range</b>	Fridge: +2° C to +10° C, Freezer: -18° C to -25° C
<b>Colour</b>	White/Green
<b>Door Type</b>	Solid
<b>Door Style</b>	Hinged, Reversible
<b>Door Open Depth (mm)</b>	1090
<b>Replacement Door Seal</b>	No
<b>Open Door Sensor</b>	Yes
<b>Lock</b>	Yes
<b>Adjustable Feet</b>	2 Adjustable Feet, 2 Rear Rollers
<b>Temperature Display on Panel</b>	Yes
<b>Temperature Alarm</b>	Yes
<b>Number of Shelves</b>	2
<b>Number of Shelf Positions</b>	4
<b>Shelf Size</b>	1: 308 x 448 x 4.5, 1: 200 x 448 x 4.5
<b>Shelf Material</b>	Wire
<b>Moulded Drawers</b>	1 Solid
<b>Adjustable Door Bins</b>	2 full + 2 half
<b>Interior Light</b>	No
<b>Fan</b>	No
<b>Spark Free</b>	Yes
<b>Auto or Manual defrost</b>	Fridge: Automatic, Freezer: Manual
<b>Controller type</b>	Digital, External
<b>Climate Class</b>	N
<b>Refrigerant Type</b>	R600a
<b>CFC Free</b>	Yes
<b>HFC Free</b>	Yes
<b>HCFC Free</b>	Yes
<b>UK 3 Pin Plug</b>	Yes
<b>Other Plug Types Available</b>	Yes
<b>Voltage Supply</b>	230
<b>Aux Output</b>	Yes
<b>Insulation Thickness (mm)</b>	45
<b>Insulation Type</b>	Cyclopentane
<b>Anti Microbial Handle</b>	Yes