

**Level 0**  
0 - <10 min

**Level 1**  
10 - 30 min

**Level 2**  
30 - 60 min

**Level 3**  
60 - 120 min

**Level 4**  
120 - 240 min

**Level 5**  
240 - 480 min

**Level 6**  
< 480 min

## duoSHIELD\* PFT Nitrile 240

CE  
0120



- Class 1 MDD glove (Council Directive 93/42/EEC)
- Category III PPE glove (PPE Directive 89/686/EEC)
- Complex Design - For mortal and irreversible risks
- Powder-free white nitrile glove
- Ambidextrous
- 240 mm / 0.08 mm (EN 455-2:2015 and EN 420:2003+A1:2009)
- Micro-organism resistant - AQL <0.65 (EN 374-2:2014 Level 3 and EN 455-1:2000)
- Virus resistant (ISO 374-5:2016 / ISO 16604:2004 Procedure B)
- Tested for chemical permeation - Type C according to ISO 374-1:2015
- EN 16523-1:2015 chemical testing method
- Type I hypersensitivity eliminated - Type IV hypersensitivity reduced

1310-73-2  
Sodium Hydroxide 40%

**LEVEL 6**  
**480 min**  
**DR 13%**

7722-84-1  
Hydrogen Peroxide 30%

**LEVEL 0**  
**6 min**  
**DR 1%**