

New and improved Portable Cooler range

With our new Portable Vaccine Coolers, you have the **versatility and simplicity** to take your vaccines with you on the move.

With two sizes ranging from **8L to 18L**, you can plan your movements for the day without the worry of cold chain compromise. Alongside this, the use of gel packs rather than a refrigeration system, means you can take this product anywhere **without the worry of needing to find a plug**.



8 Litre Cooler
VCP8UK | 444410882



18 Litre Cooler
VCP18UK | 444410883



- +2°C to +8°C operating temperature range
- External LED temperature display
- Holds temperature inside between 2°C & 8°C for up to 30 hours¹
- Comes with 9 Gel Packs (VPC8UK) and 12 Gel Packs (VPC18UK)
- Aluminium heat preserving foil lining
- Typical vaccine storage 2.5-3L (8L) and 12-13L (18L)
- Covered bung non slip feet
- Webbing adjustable shoulder strap²

1. 8L model holds the same temperature for up to 17 hours 2. 18L model only.