



Appleton PCR Plate Compatibility Chart



appleton®

For the Life Scientist

	Note/ Holder number	ACP122	ACP123	ACP024	ACP025
Standard Well (0.2ml, height 21mm.)		•			
Low Profile Well (0.1ml, height 16mm.)			•	•	•
Applied Biosystems					
7000		•			
7300	4367032	•			
7500	4367032	•			
7900HT Fast (with tray) (0.1ml type)			•		
7900HT Fast (with tray) (0.2ml type)		•			
StepOne	48-Well				
StepOne Plus			•		
7500Fast	4359652		•		
ViiA7 Fast (0.1ml model)	4459846		•		
ViiA7 Fast (0.2ml model)	4459845	•			
QuantStudio Fast 96-well Block (0.1ml model)	4459846		•		
QuantStudio Fast 96-well Block (0.2ml model)	4459845	•			
5700	Sequencer	•			
7700	Sequencer	•			
310/3100/3130	Sequencer	•			
3700/3730	Sequencer	•			
Roche					
LightCycler 480					•
LightCycler 480 II					•
LightCycler 96					•
LightCycler Nano	32-Well				
Bio-Rad					
iCycler iQ					
iQ5					
MyiQ					
Chromo4					
Opticon I				•	•
Opticon II				•	•
mini Opticon	48-Well				
CFX96				•	•
CFX384					
Stratagene					
Mx4000/3000P/3005P					
Eppendorf					
MasterCycle ep Realplex					
Thermo Scientific					
PikoReal 24 Real-Time PCR	24-Well				
Biometra					
Toptical Thermocycler					

Appleton qPCR products corresponding to Applied Biosystems & Bio-Rad Plasticware

Appleton Code	Applied Biosystems Code	Bio-Rad Code
ACP122	N8010560 / 4316813	MSS-9601
ACP123	4346907	N/A
ACP024	N/A	MSP-9601 / MLL-9601
ACP025	N/A	MLL-9651



0121 458 7740



info@appletonwoods.co.uk



appletonwoods.co.uk

appleton[®]

For the Life Scientist